

**United States Department of Agriculture
Agricultural Marketing Service, Science & Technology
Microbiological Data Program**

SOP No.: MDP-MTH-04		Page 1 of 5
Title: Detection of <i>Salmonella</i> in Fresh Produce by BAX [®] PCR		
Revision: 02	Replaces: 01	Effective: 01/01/06

1. Purpose

To provide a standard procedure for detection of *Salmonella* species in fresh produce using the BAX[®] system by all laboratories participating in the USDA/AMS Microbiological Data Program (MDP). The BAX[®] system is used here as a screening tool.

2. Scope

This standard operating procedure (SOP) shall be followed by all laboratories conducting microbiological studies for MDP, including support laboratories conducting non-routine activities. This SOP represents minimum MDP requirements and is presented as a general guideline. Each laboratory shall have written procedures that provide specific details concerning how the procedure has been implemented in that laboratory.

3. Principle

The BAX[®] PCR system is a DNA-based screening method developed by DuPont Qualicon for detecting bacterial pathogens in food and environmental samples. The sensitivity and the accuracy are a result of the use of polymerase chain reaction (PCR) to amplify DNA fragments specific for a given target organism. The amplification of specific DNA fragments is monitored using a dye that fluoresces upon intercalating with the double-stranded DNA of the amplified PCR product.

4. Outline of Procedures

Equipment and Materials	6.1
List of Controls	6.2
BAX [®] Analysis	6.3
BAX [®] -Positive Samples	6.4
Reporting	6.5

**United States Department of Agriculture
Agricultural Marketing Service, Science & Technology
Microbiological Data Program**

SOP No.: MDP-MTH-04		Page 2 of 5
Title: Detection of <i>Salmonella</i> in Fresh Produce by BAX [®] PCR		
Revision: 02	Replaces: 01	Effective: 01/01/06

5. References

- 5.1. BAX[®] System User Guide & Protocol Summary, DuPont Qualicon
- 5.2. SOP MDP-DATA-01 Record Keeping and Results Reporting
- 5.3. SOP MDP-LABOP-02, Sample Receipt, Elution, Preenrichment, and DNA Extraction
- 5.4. SOP MDP-MTH-03A, Isolation and Identification of *Salmonella* from Fresh Produce
- 5.5. SOP MDP-QA-03, Quality Assurance (QA) Controls
- 5.6. Evaluation of Enrichment Ability of Universal Pre-Enrichment Broth (UPB) for *Salmonella* ser. Typhimurium and *E. coli* O157:H7 from Produce Commodities. Final study report, Division of Consolidated Laboratory Services (DCLS), Department of General Services, Commonwealth of Virginia. October 2005.

6. Procedures

6.1. Equipment and Materials

- 6.1.1. BAX[®] System
- 6.1.2. BAX[®] PCR assay kit for *Salmonella* species
- 6.1.3. Additional materials needed to perform procedure as listed in BAX[®] System User Guide & Protocol Summary

Note: Use BAX[®] lysis buffer without protease.

6.2. List of Controls (Specific strains are listed in SOP MDP-QA-03)

- 6.2.1. DNA from Universal Preenrichment Broth (UPB) (Blank)
 - 6.2.2. DNA from negative culture control from SOP MDP-LABOP-02
 - 6.2.3. DNA from positive culture control from SOP MDP-LABOP-02
 - 6.2.4. DNA from positive produce culture control from SOP MDP-LABOP-02
 - 6.2.5. BAX[®] lysis buffer (without protease)
-

**United States Department of Agriculture
Agricultural Marketing Service, Science & Technology
Microbiological Data Program**

SOP No.: MDP-MTH-04		Page 3 of 5
Title: Detection of <i>Salmonella</i> in Fresh Produce by BAX [®] PCR		
Revision: 02	Replaces: 01	Effective: 01/01/06

6.2.6. Carry all controls through this entire procedure, including any necessary cultural confirmation. If any of the controls fail to yield a satisfactory result refer to SOP MDP-QA-03.

6.3. BAX[®] Analysis

6.3.1. Use DNA extracted from the UPB preenriched cultures. Keep the DNA under refrigeration or in a cooling block until ready to use.

6.3.2. Transfer PCR tubes to the cooling block.

6.3.3. Hydrate BAX[®] reagent pellet by adding 50 µL of BAX[®] lysis buffer (without protease) to the PCR tubes.

6.3.4. Transfer a 5-µL aliquot of extracted DNA prepared from each of the UPB-preenriched samples and controls (from SOP MDP-LABOP-02) to the PCR tubes.

Note: To minimize contamination, keep samples and controls separate.

Note: Do not add DNA sample directly to the BAX[®] reagent pellet without prior addition of lysis buffer (without protease) alone or as a DNA:lysis buffer mix.

6.3.5. Proceed as directed by BAX[®] Users manual.

6.4. BAX[®]-Positive Samples

6.4.1. When a sample is positive or indeterminate via BAX[®], proceed to attempt to isolate *Salmonella* as directed in SOP MDP-MTH-03A.

6.5. Reporting

6.5.1. A BAX[®]-positive result is considered a preliminary positive result.

6.5.2. Data shall be reported according to SOP MDP-DATA-01.

Disclaimer: Reference to brand names (kits, equipment, media, reagents, etc.) does not constitute endorsement by this agency.

**United States Department of Agriculture
Agricultural Marketing Service, Science & Technology
Microbiological Data Program**

SOP No.: MDP-MTH-04		Page 4 of 5
Title: Detection of <i>Salmonella</i> in Fresh Produce by BAX [®] PCR		
Revision: 02	Replaces: 01	Effective: 01/01/06

John Punzi *12/9/05*

Revised by: John Punzi
Chemist, Monitoring Programs Office
8609 Sudley Road, Suite 206
Manassas, VA 20110
(703) 330-2300

Date

Cindy Koschmann *12/8/05*

Approved by: Cindy Koschmann
MDP Technical Advisory Group
Wisconsin Department of Agricultural, Trade and Consumer Protection
Bureau of Lab Services
4702 University Avenue
Madison, WI 53707-7883
(608) 267-3510

Date

Diana Haynes *12/9/05*

Approved By: Diana Haynes
Deputy Director, Monitoring Programs Office
8609 Sudley Road, Suite 206
Manassas, VA 20110
(703) 330-2300

Date

*ORIGINAL SIGNATURE PAGE MAINTAINED BY USDA, AMS, SCIENCE & TECHNOLOGY, MONITORING PROGRAMS OFFICE
ELECTRONICALLY REPRODUCED SIGNATURES*

**United States Department of Agriculture
Agricultural Marketing Service, Science & Technology
Microbiological Data Program**

SOP No.: MDP-MTH-04		Page 5 of 5
Title: Detection of <i>Salmonella</i> in Fresh Produce by BAX [®] PCR		
Revision: 02	Replaces: 01	Effective: 01/01/06

Revision 02 January 2006 Monitoring Programs Office

- Added references to UPB
- Modified BAX[®] procedure to allow extracted DNA
- Screening procedure now performed on individual samples (rather than 3 pooled samples)
- All cultural procedures now included in SOP MDP-MTH-03A

Revision 01 January 2005 Monitoring Programs Office

- Removed references to outdated control strains
 - Removed all duplicative instructions found in BAX[®] Users Manual
 - Added provision to obtain 25 mL wash eluate for positive produce control
-
-