

UNITED STATES DEPARTMENT OF AGRICULTURE
AGRICULTURAL MARKETING SERVICE
FRUIT AND VEGETABLE DIVISION
PROCESSED PRODUCTS STANDARDIZATION AND INSPECTION BRANCH
WASHINGTON, D. C.

FILE CODE
135 - A - 3

**TECHNICAL
INSPECTION PROCEDURES**

BRIX MEASUREMENT

FOR USE OF USDA PROCESSED FOODS INSPECTORS

APRIL, 1960

ACTION BY: All Processed Food Inspectors

APPROVED BY: F. L. Southerland, Chief of the Branch

FILE UNDER: Chapter X - Laboratory Procedures

THESE INSTRUCTIONS SUPERSEDE:

- 1) Memorandum A-11 (11-9-49)
- 2) Memorandum A-76 (9-28-51)

TABLE OF CONTENTS

| <u>TOPIC:</u> | <u>PAGE</u> |
|--|-------------|
| I GENERAL | |
| A TERMINOLOGY | 1 |
| B FACTORS AFFECTING BRIX IN CANNED FRUITS AND SIMILAR PRODUCTS | 1 |
| II PREPARATION OF SAMPLES | |
| A EQUALIZATION - TERMINOLOGY | 1 |
| B PROCEDURE | 2 |
| 1 Equalized Products | 2 |
| 2 Non-equalized Products. | 3 |
| a Simulated Equalization. | 3 |
| b Equipment | 3 |
| c Procedure. | 4 |
| III SELECTION OF MEASURING INSTRUMENT | |
| A HYDROMETER. | 5 |
| B REFRACTOMETER. | 5 |
| IV MEASUREMENT TECHNIQUE | |
| A BRIX HYDROMETER | |
| 1 Description. | 6 |
| 2 Standardizing the Hydrometer. | 7 |
| 3 Procedure. | 8 |
| 4 Care of the Hydrometer. | 9 |
| B REFRACTOMETER | |
| 1 Description. | 10 |
| 2 Standardizing the Refractometer. | 10 |
| 3 Procedure. | 12 |
| 4 Care of the Refractometer. | 13 |
| V REPORTING RESULTS. | 14 |
| ILLUSTRATION. | 15 |
| TABLE A REFRACTIVE INDEX, DEGREES BRIX, SPECIFIC GRAVITY, AND DEGREES BAUME OF SUGAR (SUCROSE) SOLUTIONS. . . | R1-R5 |
| TABLE B TEMPERATURE CORRECTIONS - BRIX HYDROMETER. | R6,R7 |
| TABLE C TEMPERATURE CORRECTIONS - REFRACTOMETER. | R8 |

INSPECTORS' INSTRUCTIONS FOR BRIX MEASUREMENT

- APRIL 1960 -

I GENERAL

A TERMINOLOGY

The term "Brix" technically means the percent by weight of sugar solids in a pure sucrose solution. In canned fruits, fruit juice, and similar products it is commonly used as a convenient term to express the percent by weight of all soluble solids in solutions that are generally not pure sucrose but which contain sweeteners or mixtures of two or more types of sweeteners and small amounts of other substances.

The instruments generally used for Brix measurement are the Brix hydrometer and the refractometer. These instruments are affected not only by the sugars present in a product but also by such substances as fruit acids, pectins, and minerals. Regardless of the composition of the solution, the Brix reading expresses the soluble solids content of such solution in terms of a Brix value corresponding to a pure sucrose solution of the same specific gravity.

B FACTORS AFFECTING BRIX IN CANNED FRUIT AND SIMILAR PRODUCTS

The final Brix measurement of liquid packing media in canned fruits and similar products depends largely on four factors:

- (1) The soluble solids of the in-going fresh product.
- (2) The Brix of the in-going liquid medium.
- (3) The proportion of fresh product to packing medium, and
- (4) The extent to which the finished product (fresh product and liquid medium) has equalized.

II PREPARATION OF SAMPLES

A EQUALIZATION - TERMINOLOGY

In order to properly evaluate Brix readings it is important that the product and any included packing medium have undergone equalization. In the case of canned fruits

and vegetables packed in sirup, complete equalization may require several weeks. However, for the purposes of inspection, it may be assumed that equalization is completed 15 days after packing unless otherwise stated in the particular standard or specification. This, of course, does not mean that the product and packing media have completely equalized during this period but equalization is sufficiently complete for inspection and certification purposes.

It is not necessary to delay inspection pending the 15 day equalization period. Using proper equipment and techniques, complete equalization may be accomplished at any time by comminuting the product and packing medium. This procedure is referred to as "simulated equalization" and will yield comparable results to natural equalization.

In accordance with the foregoing discussion there are two basic categories of products on which Brix is measured; 1) naturally equalized products (either by prolonged storage or by the nature of the product), 2) non-equalized products.

Samples are to be taken from individual containers only. Composite samples are not to be used.

B PROCEDURE

1 Equalized Products

- a Mix the sample. Completely liquid products should be stirred while in their original container.

Liquid packing media of products such as canned fruits is recovered during or after the drain weight operation, while the fruit remains on the sieve. The grading tray, over which the individual container is drained, must be free from water and reasonably free from sirup adhering as a result of inspecting a previous sample. (Pans used for sirup samples that are borderline or in dispute shall, of course, be cleaned and dried each time). Sirup collected in this manner is considered well mixed. Liquid packing media of products such as canned fruits in

smaller containers (where the contained volume of sirup is less than or approximately equals the capacity of the glass cylinder) may be poured from the can directly into a glass cylinder when the Brix hydrometer is to be used as the measuring device.

- b Transfer sirup to cylinder and take the Brix spindle reading.
- c Make temperature corrections (if any).
- d If refractometer is used, follow procedure as outlined in 2 c, (3)-(5) of this section.

2 Non-Equalized Products

- a Simulated Equalization is the process of simulating the equalization of soluble solids in the fruit or vegetable ingredient and the packing media. For all practical purposes, this process produces the same results as those obtained by the Brix hydrometer method on the liquid packing media of the same product after complete equalization.

In essence, the entire contents of the canned products are comminuted into a homogeneous liquid state and degree Brix of this slurry is determined refractometrically. The resultant degree Brix is afforded the same significance as that obtained by Brix hydrometer on the liquid packing media of the equalized product.

b Equipment

- (1) Comminutor such as "Waring Blendor," "Osterizer," or other similar machine with bowl size adequate to accomodate contents of No. 2 1/2 can size and smaller.
- (2) Comminutor, as above, with bowl size adequate to accomodate entire contents of a No. 10 can size.
- (3) Refractometer, water cooled, and equipped with thermometer, and suitable light source.

- (4) Temperature corrections chart for refractometer Brix scale. (Included at end of this instruction).
- (5) Soft applicator such as plastic spatula or glass rod tipped with rubber policeman.
- (6) Filtering equipment:
 - (a) funnel
 - (b) test tube
 - (c) cover glass
 - (d) filter paper (such as Whatman #12).

NOTE: Nylon cloth or similar material may be used in lieu of filter paper provided comparable results can be obtained.

c Procedure

- (1) Pour entire contents of container into comminutor bowl. Whole, unpitted fruit should be hand-pitted before comminuting. Take care that all the fruit and liquid is transferred to the bowl insofar as is practicable.
- (2) Comminute the sample until homogeneous. This process should take from one to two minutes. Check efficiency of blending by occasionally pouring a sample onto a clean, dry pan and visually examining for lumps of fruit or vegetable flesh. Further mixing is indicated if these lumps are present. The inspector will arrive at correct mixing times with experience. (It is usually not possible to completely liquify maraschino cherries used in fruit cocktail and fruit mixes; however, final results are not adversely affected). Avoid excessive mixing since it shortens the motor-life of the comminutor and over heats the sample. Special comminutor blades (such as the "Cenco-Pinto" available from Central Scientific Co.) may be used where the product is packed in very small containers and difficulty is encountered with regular comminutor blades. An alternative method for overcoming any small container

problem is the use of double-dilution with an equal part by weight of distilled water, comminuting and multiplying the result by 2. Also, a small fruit jar may be adapted with blending knives and used for comminuting small samples.

- (3) Filter if necessary, discarding first few drops. Avoid evaporation.
- (4) Apply about two drops of sample to the refractometer prisms.
- (5) Take the reading, make temperature correction (if any) and record results.

III SELECTION OF MEASURING INSTRUMENT

A Hydrometer

- 1 Naturally equalized canned fruits and canned vegetables (such as sweet potatoes in sirup) consisting of separable units in a liquid packing medium, (including water and dietetic packs).
- 2 Canned, chilled and frozen single-strength fruit juices.
- 3 Frozen fruits. In-going sirups only.
- 4 Special products (such as pickles, pickle relish, molasses and table sirups) as specified in standards or specifications and instructions.

B Refractometer

- 1 Non-equalized canned fruits and canned vegetables (such as sweet potatoes in sirup) consisting of separable units in a liquid packing medium, including water and dietetic pack--after undergoing simulated equalization.
- 2 Canned fruit products such as purees, pulps, butters, jams, preserves and jellies.
- 3 Canned and frozen fruit juices other than single-strength.

- 4 Naturally equalized canned fruits and canned vegetables (such as sweet potatoes in sirup) consisting of separable units in a liquid packing medium, where individual containers contain insufficient liquid medium to measure by hydrometer.
- 5 Frozen fruits.
- 6 Special products (such as honey) as specified in standards, specifications or instructions.

IV MEASUREMENT TECHNIQUE

A BRIX HYDROMETER

1 Description

The Brix hydrometer or spindle is a glass instrument consisting of a hollow bulb which causes the hydrometer to float. Below this bulb is a section filled with shot to cause the instrument to float in an upright position, and above the bulb is a graduated stem from which the readings are taken. The extent to which the stem extends above the surface of the liquid depends upon the specific gravity of the liquid being tested.

The instrument is similar to any specific gravity hydrometer indicating degrees of liquid density. However, the graduations on the stem of this instrument are expressed in degrees Brix and indicate the percentage of soluble solids in terms of sucrose equivalent. If the Brix hydrometer is immersed in a solution composed entirely of sucrose and water, the reading at proper temperature is exactly the percentage by weight of pure dry sugar (sucrose) in the solution. The degree of Brix is easy to determine and is a convenient measurement of the specific gravity in terms of the Brix scale and reflects the total soluble solids in solution.

The hydrometer most often used in laboratories has a 5 degree range for each spindle with graduations in 1/10 degrees. Thus six or seven hydrometers are necessary to cover the range

of sirups normally encountered during inspections. Most laboratories should be equipped with instruments to cover at least the range from 5 degrees to 40 degrees Brix.

2 Standardizing the Hydrometer

Brix hydrometers should be checked periodically for accuracy. Two methods are suggested:

- a Comparison of the hydrometer reading of a sample of sirup with the reading of another hydrometer of known accuracy on the same sirup sample; or,
- b Testing the hydrometer in a specially prepared sucrose solution which contains a known percentage, by weight, of pure, dry sucrose. For all practical purposes a good grade of dry commercial cane or beet sugar will be sufficiently pure for standardization of hydrometers. Observe the following rules:
 - (1) Select a point approximately midway on the spindle to be tested. For example, a hydrometer in the range of 15 to 20 degrees can be tested in a sirup of 18°.
 - (2) Prepare a sugar solution of the required density by dissolving a weighed quantity of dry, free-running sugar in a weighed quantity of distilled water. For example, assume it is desired to prepare 1000 grams of 18° sirup. Tare a beaker or suitable receptacle on a gram balance and add exactly 180 grams of dry sugar. Then add sufficient water to bring the weight of the mixture (sugar and water) to exactly 1000 grams. This solution will give a reading of 18° Brix with an accurate hydrometer at the proper temperature. Do not make up to volume in a volumetric flask.

Bear in mind that 18° Brix is 18% by weight of sugar solids which is 18 grams per 100 grams of solution and not 18 grams made up to 100 ml.

- (3) Mix the sugar and water well. Finally transfer the contents back and forth between two large beakers in order to assure proper mix.
- (4) Cool solution down to about 1 degree below the temperature at which the instrument is calibrated.
- (5) Transfer sirup to a cylinder and observe temperature at which the instrument was calibrated. If necessary, immerse in a water bath of the proper temperature.
- (6) When the sirup in the cylinder has equalized to the proper temperature, immerse the hydrometer and observe the reading. It is well to take a series of readings by removing the hydrometer and cleaning and drying it before each immersion.
- (7) If the hydrometer reads within 1/10 of a degree of the proper reading, it can be considered sufficiently accurate for our purposes. If the error is more than 1/10 degree, a proper correction should be recorded for such error and each reading with this hydrometer so adjusted.
- (8) If the error in the instrument exceeds 3/10 degree, replace the instrument.
- (9) For more precise calibration the instrument should be standardized at more than one check point.

3 Procedure

- a Pour the liquid into a standard glass laboratory cylinder (approximately 1 1/4 inches inside diameter, eight to ten inches tall). Completely fill the cylinder, if sufficient liquid is available, and allow to overflow slightly so as to float off any foam or bubbles.

- b Slowly lower clean, dry hydrometer into the liquid until it is very near floating position. Release spindle with a slight spinning motion. If it is dropped into the sirup it will sink too far and when it rises, some of the sirup will adhere to the stem. The weight of this adhering sirup will cause the hydrometer to float at a lower position than it should. It is also important to note that air bubbles in the sirup will cause the hydrometer to float at an erroneous level.
- c Observe the reading on the stem after allowing the hydrometer to come completely to rest. The hydrometer must not touch the side of the glass cylinder.

The eye should be on a level with the surface of the liquid (see illustration). Where the liquid touches the stem, it rises a short distance to form a meniscus. The reading will be inaccurate if taken at the top of the meniscus rather than at the true liquid level. The meniscus layer is usually thin and the graduation marks can be seen through it. However, with very dark colored sirups it may be necessary to make a slight allowance for the meniscus factor.

- d Make any appropriate temperature correction. Record the temperature of the sirup and refer to the temperature corrections chart for the Brix hydrometer. The chart is self-explanatory but it should be noted that the correction factor varies according to the specific gravity of the liquid as well as with temperature. The correction factor is taken from the column nearest the degree Brix found upon examination. For example, if the hydrometer indicates a reading of 23.5 degrees Brix, the correction factor is taken from the vertical column headed "25 degrees Brix."

4 Care of the Hydrometer

The Brix hydrometer is a very delicate and fragile instrument. It must be handled with care and respect to avoid unnecessary breakage and is to be thoroughly cleaned and dried immediately after each use. Protect it from breakage when not in use by storing in an adequately cushioned container.

B REFRACTOMETER

1 Description

The refractometer is an instrument that optically measures the density of a liquid. Light is passed through the sample and is deflected in relation to the density of the sample. The instrument is calibrated in terms of refractive index and also usually contains a scale in terms of degrees Brix. A common type of refractometer consists of two prisms between which a portion of the material to be tested is placed; a mirror which reflects light through the prisms and liquid; a telescopic tube with cross-hairs superimposed on the field of vision; a scale calibrated in terms of refractive index, degrees Brix, or both; a compensator with which to correct for the chromatic dispersion of light; and a thermometer in the water circulating system.

Further descriptions of the refractometer vary with the make and model of the individual instrument. The manual of instructions accompanying each refractometer should be studied closely. Many of the small, hand-type refractometers are not sufficiently accurate for laboratory purposes but may be used to good advantage in estimating the soluble solids content of fresh fruit.

2 Standardizing the Refractometer

The refractometer should be checked for accuracy before use. It should be checked daily when in constant use, such as on In-plant inspections.

The method of determining the instrument's accuracy of calibration should be that specified in the manual of instructions accompanying the refractometer. This is generally accomplished through the use of a special prism of known refractive index.

Standardizing the refractometer by the test prism is considered the most accurate method. This method is to be used periodically.

-11-

July 1969

Frequent day-to-day checks may be made with distilled water provided there is very close correlation with the test prism method. The distilled water checks shall supplement less frequent checks with the test prism.

The distilled water check for accuracy shall be made as follows:

- (1) Adjust prism temperature as close as possible to standard temperature. (In hot weather it may be necessary to surround the water hoses with ice water, or draw from a reservoir of cooled water).
- (2) Clean and dry the prisms carefully.
- (3) Apply one or two drops of distilled water to prisms.
- (4) Compare observed refractive index and prism temperature with the appropriate refractive index (on the table below) for the same temperature. If the refractometer differs more than plus or minus 0.0002 repeat steps (2) and (3) until a total of three readings is obtained. If the refractive index of the last reading still differs from the table more than plus or minus 0.0002, check the refractometer with the standard test prism. If the reading still differs more than plus or minus 0.0002 recalibrate the refractometer.

CAUTION:

Recalibration of the refractometer should be done by only an experienced inspector.

REFRACTIVE INDEX OF DISTILLED WATER

| Temperature | | Abbe' | | Temperature | | Abbe' | |
|-------------|----|---------------|---------|-------------|----|---------------|---------|
| °F | °C | Refractometer | Reading | °F | °C | Refractometer | Reading |
| 59.0 | 15 | 1.3334 | : | 71.6 | 22 | 1.3328 | |
| 60.8 | 16 | 1.3334 | : | 73.4 | 23 | 1.3327 | |
| 62.6 | 17 | 1.3333 | : | 75.2 | 24 | 1.3326 | |
| 64.4 | 18 | 1.3332 | : | 77.0 | 25 | 1.3325 | |
| 66.2 | 19 | 1.3331 | : | 78.8 | 26 | 1.3324 | |
| 68.0 | 20 | 1.3330 | : | 82.4 | 28 | 1.3322 | |
| 69.8 | 21 | 1.3329 | : | 86.0 | 30 | 1.3320 | |

3 Procedure

- a Apply one or two drops of the sample (filtered when necessary) to the lower prism. Avoid scratching the soft, refractometer prisms.
- b Close and lock the prism chamber.
- c Apply the light source and view the shadow through the telescope, keeping the eye in the center of the eyepiece. Align the shadow edge (which should be very sharp) exactly with the intersection of the cross-hairs. An indistinct shadow edge may be caused by:
 - (1) Insufficient or excessive sample on prisms.
 - (2) Unfiltered sample of purees and slurries.
 - (3) Improper chromatic dispersion adjustment (indicated by a blue or red tinge in the shadow edge).
 - (4) Insufficient or excessive light directed on the prisms.
 - (5) Improper closing of prisms (too rapid closure of prisms or prisms not locked securely).
 - (6) Light stop (present on some instruments) that controls amount of light transmitted to prisms not fully closed.
- d Read the degree Brix from the Brix scale and note the temperature indicated by the thermometer.
- e Make any appropriate temperature correction.

Temperature greatly influences the Brix reading and therefore it is imperative that the reading be corrected to the instrument's "standard" temperature - usually 20° C (68° F). This

correction is made in the same manner as the correction adjustment for the Brix hydrometer except that a different chart is consulted. The appropriate chart is supplied at the end of this instruction. Temperature of the prisms should remain as close to standard temperature, usually 20° C (68° F), as is practical by adjustment of the flow of water through the cooling system.

CAUTION: The inclusion of water at any point in the determination will result in erroneous readings. Placing a surplus of sample on the prisms may cause a portion of the liquid to flow around onto the backside of a prism which will result in an erroneous reading.

These instructions are to be supplemented by the manual of instructions pertaining to the particular refractometer being used.

4 Care of the Refractometer

The instrument must be kept clean and in proper working order at all times.

Clean and dry the prisms with soft, lintfree tissues immediately after each reading. The prisms are very soft and scratch easily. Therefore, never touch them with hard objects such as spoons or glass rods. It is well to place a small piece of soft paper or tissue between them when the instrument is not in use.

The refractometer should be protected with a cover when not in use and should be kept in its case during long periods of storage.


Instruments standardized at 20° C are preferred by the Branch. When sample temperatures are above or below the temperature at which the instrument is calibrated, corrections are made to the standard temperature of the instrument. Temperature corrections charts included in this instruction are for instruments standardized at 20° C only.

V REPORTING RESULTS

Unless otherwise specified, report results of Brix readings to the closest 1/10th degree.

In reporting results, be sure to adjust the readings for any temperature or instrument correction.

Corrections for presence of fruit acids, minerals,
~~and similar ingredients are not to be made unless~~
~~so specified in standards, specifications or~~
Inspectors' Instructions.



-R1-

TABLE A -- REFRACTIVE INDEX, DEGREES BRIX, SPECIFIC GRAVITY
AND DEGREES BAUME' OF SUGAR (SUCROSE) SOLUTIONS

| Refractive Index at 20°C. | Degrees Brix | Specific Gravity 20°/20°C. | Degrees Baume' Mod.145 | Refractive Index at 20° C. | Degrees Brix | Specific Gravity 20°/20°C. | Degrees Baume' Mod.145 |
|---------------------------------|-----------------|----------------------------------|------------------------------|----------------------------------|-----------------|----------------------------------|------------------------------|
| 1.3330 | 0.0 | 1.00000 | 0.00 | 1.3423 | 6.4 | 1.02530 | 3.57 |
| 1.3333 | 0.2 | 1.00078 | 0.11 | 1.3426 | 6.6 | 1.02611 | 3.69 |
| 1.3336 | 0.4 | 1.00156 | 0.22 | 1.3429 | 6.8 | 1.02692 | 3.80 |
| 1.3339 | 0.6 | 1.00233 | 0.34 | 1.3432 | 7.0 | 1.02773 | 3.91 |
| 1.3341 | 0.8 | 1.00312 | 0.45 | 1.3435 | 7.2 | 1.02854 | 4.02 |
| 1.3344 | 1.0 | 1.00390 | 0.56 | 1.3438 | 7.4 | 1.02936 | 4.13 |
| 1.3347 | 1.2 | 1.00468 | 0.67 | 1.3441 | 7.6 | 1.03017 | 4.24 |
| 1.3350 | 1.4 | 1.00546 | 0.79 | 1.3444 | 7.8 | 1.03098 | 4.35 |
| 1.3353 | 1.6 | 1.00624 | 0.90 | 1.3447 | 8.0 | 1.03180 | 4.46 |
| 1.3356 | 1.8 | 1.00702 | 1.01 | 1.3450 | 8.2 | 1.03262 | 4.58 |
| 1.3359 | 2.0 | 1.00780 | 1.12 | 1.3454 | 8.4 | 1.03344 | 4.69 |
| 1.3362 | 2.2 | 1.00859 | 1.23 | 1.3457 | 8.6 | 1.03426 | 4.80 |
| 1.3365 | 2.4 | 1.00937 | 1.34 | 1.3460 | 8.8 | 1.03508 | 4.91 |
| 1.3368 | 2.6 | 1.01016 | 1.46 | 1.3463 | 9.0 | 1.03590 | 5.02 |
| 1.3370 | 2.8 | 1.01094 | 1.57 | 1.3466 | 9.2 | 1.03672 | 5.13 |
| 1.3373 | 3.0 | 1.01173 | 1.68 | 1.3469 | 9.4 | 1.03755 | 5.24 |
| 1.3376 | 3.2 | 1.01252 | 1.79 | 1.3472 | 9.6 | 1.03837 | 5.35 |
| 1.3379 | 3.4 | 1.01331 | 1.90 | 1.3475 | 9.8 | 1.03920 | 5.46 |
| 1.3382 | 3.6 | 1.01410 | 2.02 | 1.3478 | 10.0 | 1.04003 | 5.57 |
| 1.3385 | 3.8 | 1.01490 | 2.13 | 1.3481 | 10.2 | 1.04086 | 5.68 |
| 1.3388 | 4.0 | 1.01569 | 2.24 | 1.3484 | 10.4 | 1.04169 | 5.80 |
| 1.3391 | 4.2 | 1.01649 | 2.35 | 1.3487 | 10.6 | 1.04252 | 5.91 |
| 1.3394 | 4.4 | 1.01728 | 2.46 | 1.3490 | 10.8 | 1.04335 | 6.02 |
| 1.3397 | 4.6 | 1.01808 | 2.57 | 1.3493 | 11.0 | 1.04418 | 6.13 |
| 1.3400 | 4.8 | 1.01888 | 2.68 | 1.3497 | 11.2 | 1.04502 | 6.24 |
| 1.3402 | 5.0 | 1.01968 | 2.79 | 1.3500 | 11.4 | 1.04585 | 6.35 |
| 1.3405 | 5.2 | 1.02048 | 2.91 | 1.3503 | 11.6 | 1.04669 | 6.46 |
| 1.3408 | 5.4 | 1.02128 | 3.02 | 1.3506 | 11.8 | 1.04753 | 6.57 |
| 1.3411 | 5.6 | 1.02208 | 3.13 | 1.3509 | 12.0 | 1.04837 | 6.68 |
| 1.3414 | 5.8 | 1.02289 | 3.24 | 1.3512 | 12.2 | 1.04921 | 6.79 |
| 1.3417 | 6.0 | 1.02369 | 3.35 | 1.3515 | 12.4 | 1.05005 | 6.90 |
| 1.3420 | 6.2 | 1.02450 | 3.46 | 1.3518 | 12.6 | 1.05090 | 7.02 |

-R2-

TABLE A -- REFRACTIVE INDEX, DEGREES BRIX, SPECIFIC GRAVITY,
AND DEGREES BAUME' OF SUGAR (SUCROSE) SOLUTIONS

| Refractive Index at 20° C. | Degrees Brix | Specific Gravity 20°/20°C. | Degrees Baume' Mod. 145 | Refractive Index at 20°C. | Degrees Brix | Specific Gravity 20°/20°C. | Degrees Baume' Mod. 145 |
|----------------------------------|-----------------|----------------------------------|-------------------------------|---------------------------------|-----------------|----------------------------------|-------------------------------|
| 1.3522 | 12.8 | 1.05174 | 7.13 | 1.3661 | 21.4 | 1.08923 | 11.87 |
| 1.3525 | 13.0 | 1.05259 | 7.24 | 1.3665 | 21.6 | 1.09013 | 11.98 |
| 1.3528 | 13.2 | 1.05343 | 7.35 | 1.3668 | 21.8 | 1.09103 | 12.09 |
| 1.3531 | 13.4 | 1.05428 | 7.46 | 1.3672 | 22.0 | 1.09194 | 12.20 |
| 1.3534 | 13.6 | 1.05513 | 7.57 | 1.3675 | 22.2 | 1.09284 | 12.31 |
| 1.3537 | 13.8 | 1.05598 | 7.68 | 1.3678 | 22.4 | 1.09375 | 12.42 |
| 1.3541 | 14.0 | 1.05683 | 7.79 | 1.3682 | 22.6 | 1.09465 | 12.52 |
| 1.3544 | 14.2 | 1.05769 | 7.90 | 1.3685 | 22.8 | 1.09556 | 12.63 |
| 1.3547 | 14.4 | 1.05854 | 8.01 | 1.3688 | 23.0 | 1.09647 | 12.74 |
| 1.3550 | 14.6 | 1.05940 | 8.12 | 1.3692 | 23.2 | 1.09738 | 12.85 |
| 1.3553 | 14.8 | 1.06025 | 8.23 | 1.3695 | 23.4 | 1.09829 | 12.96 |
| 1.3556 | 15.0 | 1.06111 | 8.34 | 1.3699 | 23.6 | 1.09921 | 13.07 |
| 1.3560 | 15.2 | 1.06197 | 8.45 | 1.3702 | 23.8 | 1.10012 | 13.18 |
| 1.3563 | 15.4 | 1.06283 | 8.56 | 1.3706 | 24.0 | 1.10104 | 13.29 |
| 1.3566 | 15.6 | 1.06369 | 8.67 | 1.3709 | 24.2 | 1.10195 | 13.40 |
| 1.3569 | 15.8 | 1.06455 | 8.78 | 1.3712 | 24.4 | 1.10287 | 13.51 |
| 1.3573 | 16.0 | 1.06542 | 8.89 | 1.3716 | 24.6 | 1.10370 | 13.62 |
| 1.3576 | 16.2 | 1.06629 | 9.00 | 1.3719 | 24.8 | 1.10471 | 13.73 |
| 1.3579 | 16.4 | 1.06715 | 9.11 | 1.3723 | 25.0 | 1.10564 | 13.84 |
| 1.3582 | 16.6 | 1.06802 | 9.22 | 1.3726 | 25.2 | 1.10656 | 13.95 |
| 1.3586 | 16.8 | 1.06889 | 9.33 | 1.3729 | 25.4 | 1.10748 | 14.06 |
| 1.3589 | 17.0 | 1.06976 | 9.45 | 1.3733 | 25.6 | 1.10841 | 14.17 |
| 1.3592 | 17.2 | 1.07063 | 9.56 | 1.3736 | 25.8 | 1.10934 | 14.28 |
| 1.3595 | 17.4 | 1.07151 | 9.67 | 1.3740 | 26.0 | 1.11027 | 14.39 |
| 1.3598 | 17.6 | 1.07238 | 9.78 | 1.3743 | 26.2 | 1.11120 | 14.49 |
| 1.3602 | 17.8 | 1.07325 | 9.89 | 1.3747 | 26.4 | 1.11213 | 14.60 |
| 1.3605 | 18.0 | 1.07413 | 10.00 | 1.3750 | 26.6 | 1.11306 | 14.71 |
| 1.3608 | 18.2 | 1.07501 | 10.11 | 1.3753 | 26.8 | 1.11400 | 14.82 |
| 1.3612 | 18.4 | 1.07589 | 10.22 | 1.3757 | 27.0 | 1.11493 | 14.93 |
| 1.3615 | 18.6 | 1.07677 | 10.33 | 1.3761 | 27.2 | 1.11587 | 15.04 |
| 1.3618 | 18.8 | 1.07765 | 10.44 | 1.3764 | 27.4 | 1.11681 | 15.15 |
| 1.3621 | 19.0 | 1.07853 | 10.55 | 1.3768 | 27.6 | 1.11775 | 15.26 |
| 1.3625 | 19.2 | 1.07942 | 10.66 | 1.3771 | 27.8 | 1.11869 | 15.37 |
| 1.3628 | 19.4 | 1.08030 | 10.77 | 1.3775 | 28.0 | 1.11963 | 15.48 |
| 1.3631 | 19.6 | 1.08119 | 10.88 | 1.3778 | 28.2 | 1.12058 | 15.59 |
| 1.3635 | 19.8 | 1.08208 | 10.99 | 1.3782 | 28.4 | 1.12152 | 15.69 |
| 1.3638 | 20.0 | 1.08297 | 11.10 | 1.3785 | 28.6 | 1.12247 | 15.80 |
| 1.3641 | 20.2 | 1.08386 | 11.21 | 1.3789 | 28.8 | 1.12342 | 15.91 |
| 1.3645 | 20.4 | 1.08475 | 11.32 | 1.3792 | 29.0 | 1.12436 | 16.02 |
| 1.3648 | 20.6 | 1.08565 | 11.43 | 1.3796 | 29.2 | 1.12532 | 16.13 |
| 1.3651 | 20.8 | 1.08654 | 11.54 | 1.3800 | 29.4 | 1.12627 | 16.24 |
| 1.3655 | 21.0 | 1.08744 | 11.65 | 1.3803 | 29.6 | 1.12722 | 16.35 |
| 1.3658 | 21.2 | 1.08834 | 11.76 | 1.3807 | 29.8 | 1.12817 | 16.46 |

-R3-

TABLE A -- REFRACTIVE INDEX, DEGREES BRIX, SPECIFIC GRAVITY,
AND DEGREES BAUME OF SUGAR (SUCROSE) SOLUTIONS

| Refractive Index at 20°C. | Degrees Brix | Specific Gravity 20°/20°C. | Degrees Baume' Mod. 145 | Refractive Index at 20° C. | Degrees Brix | Specific Gravity 20°/20°C | Degrees Baume' Mod. 145 |
|---------------------------------|-----------------|----------------------------------|-------------------------------|----------------------------------|-----------------|---------------------------------|-------------------------------|
| 1.3810 | 30.0 | 1.12913 | 16.57 | 1.3970 | 38.6 | 1.17158 | 21.21 |
| 1.3814 | 30.2 | 1.13009 | 16.67 | 1.3974 | 38.8 | 1.17260 | 21.32 |
| 1.3818 | 30.4 | 1.13105 | 16.78 | 1.3978 | 39.0 | 1.17362 | 21.43 |
| 1.3821 | 30.6 | 1.13201 | 16.89 | 1.3982 | 39.2 | 1.17464 | 21.54 |
| 1.3825 | 30.8 | 1.13297 | 17.00 | 1.3986 | 39.4 | 1.17566 | 21.64 |
| 1.3829 | 31.0 | 1.13394 | 17.11 | 1.3989 | 39.6 | 1.17669 | 21.75 |
| 1.3832 | 31.2 | 1.13490 | 17.22 | 1.3993 | 39.8 | 1.17772 | 21.86 |
| 1.3836 | 31.4 | 1.13587 | 17.33 | 1.3997 | 40.0 | 1.17874 | 21.97 |
| 1.3839 | 31.6 | 1.13683 | 17.43 | 1.4001 | 40.2 | 1.17977 | 22.07 |
| 1.3843 | 31.8 | 1.13780 | 17.54 | 1.4005 | 40.4 | 1.18080 | 22.18 |
| 1.3847 | 32.0 | 1.13877 | 17.65 | 1.4008 | 40.6 | 1.18183 | 22.29 |
| 1.3850 | 32.2 | 1.13974 | 17.76 | 1.4012 | 40.8 | 1.18287 | 22.39 |
| 1.3854 | 32.4 | 1.14072 | 17.87 | 1.4016 | 41.0 | 1.18390 | 22.50 |
| 1.3858 | 32.6 | 1.14169 | 17.98 | 1.4020 | 41.2 | 1.18494 | 22.61 |
| 1.3861 | 32.8 | 1.14267 | 18.08 | 1.4024 | 41.4 | 1.18598 | 22.72 |
| 1.3865 | 33.0 | 1.14364 | 18.19 | 1.4028 | 41.6 | 1.18702 | 22.82 |
| 1.3869 | 33.2 | 1.14462 | 18.30 | 1.4032 | 41.8 | 1.18806 | 22.93 |
| 1.3872 | 33.4 | 1.14560 | 18.41 | 1.4036 | 42.0 | 1.18910 | 23.04 |
| 1.3876 | 33.6 | 1.14658 | 18.52 | 1.4040 | 42.2 | 1.19014 | 23.14 |
| 1.3879 | 33.8 | 1.14757 | 18.63 | 1.4044 | 42.4 | 1.19119 | 23.25 |
| 1.3883 | 34.0 | 1.14855 | 18.73 | 1.4048 | 42.6 | 1.19224 | 23.36 |
| 1.3887 | 34.2 | 1.14954 | 18.84 | 1.4052 | 42.8 | 1.19329 | 23.46 |
| 1.3891 | 34.4 | 1.15052 | 18.95 | 1.4056 | 43.0 | 1.19434 | 23.57 |
| 1.3894 | 34.6 | 1.15151 | 19.06 | 1.4060 | 43.2 | 1.19539 | 23.68 |
| 1.3898 | 34.8 | 1.15250 | 19.17 | 1.4064 | 43.4 | 1.19644 | 23.78 |
| 1.3902 | 35.0 | 1.15350 | 19.28 | 1.4068 | 43.6 | 1.19749 | 23.89 |
| 1.3906 | 35.2 | 1.15449 | 19.38 | 1.4072 | 43.8 | 1.19855 | 24.00 |
| 1.3909 | 35.4 | 1.15548 | 19.49 | 1.4076 | 44.0 | 1.19961 | 24.10 |
| 1.3913 | 35.6 | 1.15648 | 19.60 | 1.4080 | 44.2 | 1.20066 | 24.21 |
| 1.3916 | 35.8 | 1.15747 | 19.71 | 1.4084 | 44.4 | 1.20172 | 24.32 |
| 1.3920 | 36.0 | 1.15847 | 19.81 | 1.4088 | 44.6 | 1.20279 | 24.42 |
| 1.3924 | 36.2 | 1.15947 | 19.92 | 1.4092 | 44.8 | 1.20385 | 24.53 |
| 1.3928 | 36.4 | 1.16047 | 20.03 | 1.4096 | 45.0 | 1.20491 | 24.63 |
| 1.3931 | 36.6 | 1.16148 | 20.14 | 1.4100 | 45.2 | 1.20598 | 24.74 |
| 1.3935 | 36.8 | 1.16248 | 20.25 | 1.4104 | 45.4 | 1.20705 | 24.85 |
| 1.3939 | 37.0 | 1.16349 | 20.35 | 1.4109 | 45.6 | 1.20812 | 24.95 |
| 1.3943 | 37.2 | 1.16449 | 20.46 | 1.4113 | 45.8 | 1.20919 | 25.06 |
| 1.3947 | 37.4 | 1.16550 | 20.57 | 1.4117 | 46.0 | 1.21026 | 25.17 |
| 1.3950 | 37.6 | 1.16652 | 20.68 | 1.4121 | 46.2 | 1.21133 | 25.27 |
| 1.3954 | 37.8 | 1.16752 | 20.78 | 1.4125 | 46.4 | 1.21241 | 25.38 |
| 1.3958 | 38.0 | 1.16853 | 20.89 | 1.4129 | 46.6 | 1.21349 | 25.48 |
| 1.3962 | 38.2 | 1.16955 | 21.00 | 1.4133 | 46.8 | 1.21456 | 25.59 |
| 1.3966 | 38.4 | 1.17056 | 21.11 | 1.4137 | 47.0 | 1.21564 | 25.70 |

-R4-

TABLE A -- REFRACTIVE INDEX, DEGREES BRIX, SPECIFIC GRAVITY,
AND DEGREES BAUME' OF SUGAR (SUCROSE) SOLUTIONS

| Refractive Index at 20° C. | Degrees Brix | Specific Gravity 20°/20°C. | Degrees Baume' Mod.145 | Refractive Index at 20°C. | Degrees Brix | Specific Gravity 20°/20°C | Degrees Baume' Mod. 145 |
|----------------------------------|-----------------|----------------------------------|------------------------------|---------------------------------|-----------------|---------------------------------|-------------------------------|
| 1.4141 | 47.2 | 1.21673 | 25.80 | 1.4348 | 56.8 | 1.27041 | 30.83 |
| 1.4145 | 47.4 | 1.21781 | 25.91 | 1.4352 | 57.0 | 1.27156 | 30.94 |
| 1.4150 | 47.6 | 1.21889 | 26.01 | 1.4356 | 57.2 | 1.27272 | 31.04 |
| 1.4154 | 47.8 | 1.21998 | 26.12 | 1.4361 | 57.4 | 1.27388 | 31.15 |
| 1.4158 | 48.0 | 1.22106 | 26.23 | 1.4365 | 57.6 | 1.27504 | 31.25 |
| 1.4162 | 48.2 | 1.22215 | 26.33 | 1.4370 | 57.8 | 1.27620 | 31.35 |
| 1.4166 | 48.4 | 1.22324 | 26.44 | 1.4374 | 58.0 | 1.27736 | 31.46 |
| 1.4171 | 48.6 | 1.22434 | 26.54 | 1.4379 | 58.2 | 1.27853 | 31.56 |
| 1.4175 | 48.8 | 1.22543 | 26.65 | 1.4383 | 58.4 | 1.27969 | 31.66 |
| 1.4179 | 49.0 | 1.22652 | 26.75 | 1.4388 | 58.6 | 1.28086 | 31.76 |
| 1.4183 | 49.2 | 1.22762 | 26.86 | 1.4392 | 58.8 | 1.28203 | 31.87 |
| 1.4187 | 49.4 | 1.22872 | 26.96 | 1.4397 | 59.0 | 1.28320 | 31.97 |
| 1.4192 | 49.6 | 1.22982 | 27.07 | 1.4401 | 59.2 | 1.28437 | 32.07 |
| 1.4196 | 49.8 | 1.23092 | 27.18 | 1.4406 | 59.4 | 1.28555 | 32.18 |
| 1.4201 | 50.0 | 1.23202 | 27.28 | 1.4410 | 59.6 | 1.28672 | 32.28 |
| 1.4205 | 50.2 | 1.23313 | 27.39 | 1.4415 | 59.8 | 1.28789 | 32.38 |
| 1.4209 | 50.4 | 1.23423 | 27.49 | 1.4419 | 60.0 | 1.28908 | 32.49 |
| 1.4214 | 50.6 | 1.23534 | 27.60 | 1.4424 | 60.2 | 1.29025 | 32.59 |
| 1.4218 | 50.8 | 1.23645 | 27.70 | 1.4428 | 60.4 | 1.29143 | 32.69 |
| 1.4222 | 51.0 | 1.23756 | 27.81 | 1.4433 | 60.6 | 1.29262 | 32.79 |
| 1.4226 | 51.2 | 1.23867 | 27.91 | 1.4437 | 60.8 | 1.29380 | 32.90 |
| 1.4230 | 51.4 | 1.23978 | 28.02 | 1.4442 | 61.0 | 1.29498 | 33.00 |
| 1.4235 | 51.6 | 1.24089 | 28.12 | 1.4447 | 61.2 | 1.29618 | 33.10 |
| 1.4239 | 51.8 | 1.24201 | 28.23 | 1.4451 | 61.4 | 1.29736 | 33.20 |
| 1.4243 | 52.0 | 1.24313 | 28.33 | 1.4456 | 61.6 | 1.29855 | 33.31 |
| 1.4248 | 52.2 | 1.24425 | 28.44 | 1.4460 | 61.8 | 1.29975 | 33.41 |
| 1.4252 | 52.4 | 1.24537 | 28.54 | 1.4465 | 62.0 | 1.30093 | 33.51 |
| 1.4256 | 52.6 | 1.24649 | 28.65 | 1.4470 | 62.2 | 1.30212 | 33.61 |
| 1.4260 | 52.8 | 1.24761 | 28.75 | 1.4474 | 62.4 | 1.30334 | 33.72 |
| 1.4265 | 53.0 | 1.24874 | 28.86 | 1.4479 | 62.6 | 1.30453 | 33.82 |
| 1.4269 | 53.2 | 1.24987 | 28.96 | 1.4483 | 62.8 | 1.30537 | 33.92 |
| 1.4273 | 53.4 | 1.25099 | 29.06 | 1.4488 | 63.0 | 1.30694 | 34.02 |
| 1.4278 | 53.6 | 1.25212 | 29.17 | 1.4493 | 63.2 | 1.30815 | 34.12 |
| 1.4282 | 53.8 | 1.25325 | 29.27 | 1.4497 | 63.4 | 1.30936 | 34.23 |
| 1.4286 | 54.0 | 1.25439 | 29.38 | 1.4502 | 63.6 | 1.31055 | 34.33 |
| 1.4291 | 54.2 | 1.25552 | 29.48 | 1.4507 | 63.8 | 1.31177 | 34.43 |
| 1.4295 | 54.4 | 1.25666 | 29.59 | 1.4511 | 64.0 | 1.31297 | 34.53 |
| 1.4299 | 54.6 | 1.25780 | 29.69 | 1.4516 | 64.2 | 1.31418 | 34.63 |
| 1.4304 | 54.8 | 1.25893 | 29.80 | 1.4521 | 64.4 | 1.31540 | 34.74 |
| 1.4308 | 55.0 | 1.26007 | 29.90 | 1.4525 | 64.6 | 1.31661 | 34.84 |
| 1.4312 | 55.2 | 1.26122 | 30.00 | 1.4530 | 64.8 | 1.31784 | 34.94 |
| 1.4317 | 55.4 | 1.26236 | 30.11 | 1.4535 | 65.0 | 1.31905 | 35.04 |
| 1.4321 | 55.6 | 1.26350 | 30.21 | 1.4539 | 65.2 | 1.32028 | 35.14 |
| 1.4326 | 55.8 | 1.26465 | 30.32 | 1.4544 | 65.4 | 1.32150 | 35.24 |
| 1.4330 | 56.0 | 1.26580 | 30.42 | 1.4549 | 65.6 | 1.32271 | 35.34 |
| 1.4334 | 56.2 | 1.26695 | 30.52 | 1.4553 | 65.8 | 1.32393 | 35.45 |
| 1.4339 | 56.4 | 1.26810 | 30.63 | 1.4558 | 66.0 | 1.32516 | 35.55 |
| 1.4343 | 56.6 | 1.26925 | 30.73 | 1.4563 | 66.2 | 1.32638 | 35.65 |

-R5-

TABLE A -- REFRACTIVE INDEX, DEGREES BRIX, SPECIFIC GRAVITY,
AND DEGREES BAUME' OF SUGAR (SUCROSE) SOLUTIONS

| Refractive Index at 20°C | Degrees Brix | Specific Gravity 20°/20°C | Degrees Baume' Mod.145 | Refractive Index at 20°C. | Degrees Brix | Specific Gravity 20°/20°C | Degrees Baume' Mod.145 |
|--------------------------------|-----------------|---------------------------------|------------------------------|---------------------------------|-----------------|---------------------------------|------------------------------|
| 1.4568 | 66.4 | 1.32759 | 35.75 | 1.4798 | 75.8 | 1.38705 | 40.43 |
| 1.4572 | 66.6 | 1.32884 | 35.85 | 1.4803 | 76.0 | 1.38835 | 40.53 |
| 1.4577 | 66.8 | 1.33007 | 35.95 | 1.4808 | 76.2 | 1.38967 | 40.62 |
| 1.4582 | 67.0 | 1.33129 | 36.05 | 1.4814 | 76.4 | 1.39097 | 40.72 |
| 1.4587 | 67.2 | 1.33254 | 36.15 | 1.4819 | 76.6 | 1.39228 | 40.82 |
| 1.4591 | 67.4 | 1.33377 | 36.25 | 1.4824 | 76.8 | 1.39358 | 40.92 |
| 1.4596 | 67.6 | 1.33500 | 36.35 | 1.4829 | 77.0 | 1.39489 | 41.01 |
| 1.4601 | 67.8 | 1.33625 | 36.45 | 1.4834 | 77.2 | 1.39619 | 41.11 |
| 1.4606 | 68.0 | 1.33748 | 36.55 | 1.4839 | 77.4 | 1.39750 | 41.21 |
| 1.4611 | 68.2 | 1.33872 | 36.66 | 1.4844 | 77.6 | 1.39882 | 41.31 |
| 1.4615 | 68.4 | 1.33997 | 36.76 | 1.4849 | 77.8 | 1.40014 | 41.40 |
| 1.4620 | 68.6 | 1.34121 | 36.86 | 1.4854 | 78.0 | 1.40146 | 41.50 |
| 1.4625 | 68.8 | 1.34245 | 36.96 | 1.4860 | 78.2 | 1.40277 | 41.60 |
| 1.4630 | 69.0 | 1.34371 | 37.06 | 1.4865 | 78.4 | 1.40409 | 41.70 |
| 1.4635 | 69.2 | 1.34495 | 37.16 | 1.4870 | 78.6 | 1.40541 | 41.79 |
| 1.4640 | 69.4 | 1.34621 | 37.26 | 1.4875 | 78.8 | 1.40674 | 41.89 |
| 1.4644 | 69.6 | 1.34746 | 37.36 | 1.4880 | 79.0 | 1.40806 | 41.99 |
| 1.4649 | 69.8 | 1.34871 | 37.46 | 1.4886 | 79.2 | 1.40938 | 42.08 |
| 1.4654 | 70.0 | 1.34997 | 37.56 | 1.4891 | 79.4 | 1.41072 | 42.18 |
| 1.4659 | 70.2 | 1.35123 | 37.66 | 1.4896 | 79.6 | 1.41204 | 42.28 |
| 1.4664 | 70.4 | 1.35248 | 37.76 | 1.4901 | 79.8 | 1.41337 | 42.37 |
| 1.4669 | 70.6 | 1.35375 | 37.86 | 1.4906 | 80.0 | 1.41471 | 42.47 |
| 1.4674 | 70.8 | 1.35501 | 37.96 | | | | |
| 1.4679 | 71.0 | 1.35627 | 38.06 | | | | |
| 1.4684 | 71.2 | 1.35754 | 38.16 | | | | |
| 1.4688 | 71.4 | 1.35881 | 38.26 | | | | |
| 1.4693 | 71.6 | 1.36008 | 38.35 | | | | |
| 1.4698 | 71.8 | 1.36135 | 38.45 | | | | |
| 1.4703 | 72.0 | 1.36261 | 38.55 | | | | |
| 1.4708 | 72.2 | 1.36389 | 38.65 | | | | |
| 1.4713 | 72.4 | 1.36516 | 38.75 | | | | |
| 1.4718 | 72.6 | 1.36643 | 38.85 | | | | |
| 1.4723 | 72.8 | 1.36771 | 38.95 | | | | |
| 1.4728 | 73.0 | 1.36900 | 39.05 | | | | |
| 1.4733 | 73.2 | 1.37028 | 39.15 | | | | |
| 1.4738 | 73.4 | 1.37156 | 39.25 | | | | |
| 1.4743 | 73.6 | 1.37283 | 39.35 | | | | |
| 1.4748 | 73.8 | 1.37411 | 39.44 | | | | |
| 1.4753 | 74.0 | 1.37541 | 39.54 | | | | |
| 1.4758 | 74.2 | 1.37669 | 39.64 | | | | |
| 1.4763 | 74.4 | 1.37798 | 39.74 | | | | |
| 1.4768 | 74.6 | 1.37928 | 39.84 | | | | |
| 1.4773 | 74.8 | 1.38057 | 39.94 | | | | |
| 1.4778 | 75.0 | 1.38187 | 40.03 | | | | |
| 1.4783 | 75.2 | 1.38316 | 40.13 | | | | |
| 1.4788 | 75.4 | 1.38445 | 40.23 | | | | |
| 1.4793 | 75.6 | 1.38575 | 40.33 | | | | |

TABLE B -- TEMPERATURE CORRECTIONS -- BRUX HYDROMETER

(Standardized at 20° C.)

Based on Table 110 Bureau of Standards Circular C-440 - Nearest 1/10th Degree

Degrees Brix

Temp.
C.

| | 0 | 5 | 10 | 15 | 20 | 25 | 30 | 35 | 40 | 45 | 50 | 55 | 60 | 65 | 70 |
|----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| 0 | 0.3 | 0.5 | 0.7 | 0.8 | 0.9 | 1.0 | 1.1 | 1.2 | 1.2 | 1.3 | 1.4 | 1.4 | 1.4 | 1.5 | 1.5 |
| 5 | .4 | .5 | .6 | .7 | .7 | .8 | .9 | .9 | 1.0 | 1.0 | 1.1 | 1.1 | 1.1 | 1.1 | 1.1 |
| 10 | .3 | .4 | .4 | .5 | .5 | .6 | .6 | .6 | .7 | .7 | .7 | .7 | .8 | .8 | .8 |
| 11 | .3 | .4 | .4 | .4 | .5 | .5 | .5 | .5 | .6 | .6 | .6 | .6 | .7 | .7 | .7 |
| 12 | .3 | .3 | .4 | .4 | .4 | .4 | .4 | .4 | .5 | .5 | .5 | .5 | .6 | .6 | .6 |
| 13 | .3 | .3 | .3 | .4 | .4 | .4 | .4 | .4 | .5 | .5 | .5 | .5 | .5 | .5 | .5 |
| 14 | .2 | .3 | .3 | .3 | .3 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 |
| 15 | .2 | .2 | .2 | .3 | .3 | .3 | .3 | .3 | .3 | .4 | .4 | .4 | .4 | .4 | .4 |

-- Subtract From Observed Reading --

| | 0 | 5 | 10 | 15 | 20 | 25 | 30 | 35 | 40 | 45 | 50 | 55 | 60 | 65 | 70 |
|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|
| 16 | .2 | .2 | .2 | .2 | .2 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 |
| 17 | .1 | .1 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 |
| 18 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 |
| 19 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 |

Bo

| | 0 | 5 | 10 | 15 | 20 | 25 | 30 | 35 | 40 | 45 | 50 | 55 | 60 | 65 | 70 |
|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|----|
| 16 | .2 | .2 | .2 | .2 | .2 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 |
| 17 | .1 | .1 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 |
| 18 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 |
| 19 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 |

-- Add to Observed Reading --

| | 0 | 5 | 10 | 15 | 20 | 25 | 30 | 35 | 40 | 45 | 50 | 55 | 60 | 65 | 70 |
|----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| 21 | .0 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 |
| 22 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 |
| 23 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 |
| 24 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 |
| 25 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 |
| 26 | .3 | .3 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 |
| 27 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 |
| 28 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 |
| 29 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 |
| 30 | .6 | .6 | .6 | .6 | .6 | .6 | .6 | .6 | .6 | .6 | .6 | .6 | .6 | .6 | .6 |
| 35 | 1.0 | 1.0 | 1.0 | 1.1 | 1.1 | 1.1 | 1.1 | 1.1 | 1.1 | 1.1 | 1.1 | 1.1 | 1.1 | 1.1 | 1.1 |
| 40 | 1.4 | 1.5 | 1.5 | 1.5 | 1.5 | 1.6 | 1.6 | 1.6 | 1.6 | 1.7 | 1.7 | 1.7 | 1.7 | 1.7 | 1.7 |

TABLE B -- TEMPERATURE CORRECTIONS -- BRUX HYDROMETER (CONTINUED)

(Standardized at 20° C.)

Based on Table 110 Bureau of Standards Circular C-440 - Nearest 1/10 Degree

| | | Degrees Brix | | | | | | | | | | | | | | |
|-------------------------------|-----|--------------|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| Temp. | | 0 | 5 | 10 | 15 | 20 | 25 | 30 | 35 | 40 | 45 | 50 | 55 | 60 | 65 | 70 |
| C. | F. | | | | | | | | | | | | | | | |
| 45 | 113 | 1.9 | 1.9 | 2.0 | 2.0 | 2.0 | 2.1 | 2.1 | 2.1 | 2.1 | 2.1 | 2.1 | 2.1 | 2.1 | 2.1 | 2.1 |
| 50 | 122 | 2.5 | 2.5 | 2.5 | 2.5 | 2.6 | 2.6 | 2.6 | 2.6 | 2.6 | 2.6 | 2.6 | 2.6 | 2.6 | 2.5 | 2.5 |
| 55 | 131 | 3.1 | 3.1 | 3.1 | 3.1 | 3.1 | 3.1 | 3.1 | 3.1 | 3.1 | 3.1 | 3.1 | 3.1 | 3.0 | 3.0 | 3.0 |
| -- Add to Observed Reading -- | | | | | | | | | | | | | | | | |

-- Add to Observed Reading --

$$C = \frac{5}{9} (F - 32)$$

$$F = \left(\frac{9}{5} C \right) + 32$$

TABLE C --- TEMPERATURE CORRECTIONS -- REFRACTOMETER
(For Refractometers Standardized at 20° C.)

| Temp. ° C. | Degrees Brix | | | | | | | | | | | | | | |
|---------------|---------------------------------------|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| | 0 | 5 | 10 | 15 | 20 | 25 | 30 | 35 | 40 | 45 | 50 | 55 | 60 | 65 | 70 |
| | -- Subtract from Degrees Brix Reading | | | | | | | | | | | | | | |
| 10 | 0.5 | 0.5 | 0.6 | 0.6 | 0.6 | 0.7 | 0.7 | 0.7 | 0.7 | 0.7 | 0.7 | 0.8 | 0.8 | 0.8 | 0.8 |
| 11 | .5 | .5 | .5 | .6 | .6 | .6 | .6 | .6 | .7 | .7 | .7 | .7 | .7 | .7 | .7 |
| 12 | .4 | .5 | .5 | .5 | .5 | .5 | .6 | .6 | .6 | .6 | .6 | .6 | .6 | .6 | .6 |
| 13 | .4 | .4 | .4 | .4 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .6 | .6 |
| 14 | .3 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .5 | .5 | .5 | .5 | .5 | .5 |
| 15 | .3 | .3 | .3 | .3 | .3 | .3 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 |
| 16 | .2 | .2 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 |
| 17 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 |
| 18 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 |
| 19 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 |
| | -- Add to Degrees Brix Reading -- | | | | | | | | | | | | | | |
| 21 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 | .1 |
| 22 | .1 | .1 | .1 | .1 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 |
| 23 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 | .2 |
| 24 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 | .3 |
| 25 | .3 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 | .4 |
| 26 | .4 | .4 | .4 | .4 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 | .5 |
| 27 | .5 | .5 | .5 | .5 | .5 | .6 | .6 | .6 | .6 | .6 | .6 | .6 | .6 | .6 | .6 |
| 28 | .6 | .6 | .6 | .6 | .6 | .6 | .6 | .6 | .6 | .6 | .6 | .6 | .6 | .6 | .6 |
| 29 | .6 | .7 | .7 | .7 | .7 | .7 | .7 | .7 | .7 | .7 | .7 | .7 | .7 | .7 | .7 |
| 30 | .7 | .7 | .8 | .8 | .8 | .8 | .8 | .8 | .8 | .8 | .8 | .8 | .8 | .8 | .8 |

International Temperature Correction Table, 1936. Nearest 1/10th Degree.